

## MP68 Polyclonal Antibody

Cat No: HR1AP9170

For research use only

### Overview

Product Name	MP68 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	MP68
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MP68; C14orf2; 6.8 kDa mitochondrial proteolipid
Background	similarity:Belongs to the small mitochondrial proteolipid family.,