

MIP-T3 Polyclonal Antibody

Cat No: HR1AP6578

For research use only

Overview

Product Name	MIP-T3 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	78kDa
GeneID?Human?	TRAF3IP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TRAF3IP1; MIPT3; TRAF3-interacting protein 1; Interleukin-13 receptor alpha 1-binding protein 1; Microtubule-interacting protein associated with TRAF3; MIP-T3
Background	function:Play an inhibitory role on IL13 signaling by binding to IL13RA1. Involved in suppression of IL13-induced STAT6 phosphorylation, transcriptional activity and DNA-binding. Recruits TRAF3 and DISC1 to the microtubules.,similarity:Belongs to the TRAF3IP1 family.,subcellular location:Microtubules.,subunit:Interacts with IL13RA1. Binds to microtubules, TRAF3 and DISC1.,tissue specificity:Ubiquitous.,