

MLF1 Polyclonal Antibody

Cat No: HR1AP7358

For research use only

Overview

Product Name	MLF1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	MLF1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MLF1; Myeloid leukemia factor 1; Myelodysplasia-myeloid leukemia factor 1
Background	myeloid leukemia factor 1(MLF1) Homo sapiens This gene encodes an oncoprotein which is thought to play a role in the phenotypic determination of hemopoetic cells. Translocations between this gene and nucleophosmin have been associated with myelodysplastic syndrome and acute myeloid leukemia. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2010],