

# Megsin Polyclonal Antibody

Cat No: HR1AP8975

For research use only

## Overview

Product Name	Megsin Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	43kDa
GeneID?Human?	SERPINB7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SERPINB7; Serpin B7; Megsin; TP55
Background	serpin family B member 7(SERPINB7) Homo sapiens This gene encodes a member of a family of proteins which function as protease inhibitors. Expression of this gene is upregulated in IgA nephropathy and mutations have been found to cause palmoplantar keratoderma, Nagashima type. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2014],