

MGST3 Polyclonal Antibody

Cat No: HR1AP7352

For research use only

Overview

Product Name	MGST3 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	MGST3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MGST3; Microsomal glutathione S-transferase 3; Microsomal GST-3; Microsomal GST-III
Background	<p>microsomal glutathione S-transferase 3(MGST3) Homo sapiens This gene encodes a member of the MAPEG (Membrane Associated Proteins in Eicosanoid and Glutathione metabolism) protein family. Members of this family are involved in the production of leukotrienes and prostaglandin E, important mediators of inflammation. This gene encodes an enzyme which catalyzes the conjugation of leukotriene A4 and reduced glutathione to produce leukotriene C4. This enzyme also demonstrates glutathione-dependent peroxidase activity towards lipid hydroperoxides.[provided by RefSeq, May 2011],</p>