

IFIT-5 Polyclonal Antibody

Cat No: HR1AP6504

For research use only

Overview

Product Name	IFIT-5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GeneID?Human?	IFIT5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	IFIT5; ISG58; RI58; Interferon-induced protein with tetratricopeptide repeats 5; IFIT-5; Interferon-induced 58 kDa protein; Retinoic acid- and interferon-inducible 58 kDa protein; P58
Background	induction:By interferons.,similarity:Belongs to the IFIT family.,similarity:Contains 8 TPR repeats.,