

GCNT7 Polyclonal Antibody

Cat No: HR1AP5936

For research use only

Overview

Product Name	GCNT7 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	49kDa
GeneID?Human?	GCNT7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GCNT7; C20orf105; Beta-1; 3-galactosyl-O-glycosyl-glycoprotein beta-1,6-N-acetylglucosaminyltransferase 7
Background	function:Glycosyltransferase.,pathway:Protein modification; protein glycosylation.,sequence caution:Translated as Gln.,similarity:Belongs to the glycosyltransferase 14 family.,