

Galactosidase ?1L3 Polyclonal Antibody

Cat No: HR1AP5744

For research use only

Overview

Product Name	Galactosidase ?1L3 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	GLB1L3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GLB1L3; Beta-galactosidase-1-like protein 3
Background	similarity:Belongs to the glycosyl hydrolase 35 family.,