

## Gab 2/4 Polyclonal Antibody

Cat No: HR1AP5854

For research use only

### Overview

Product Name	Gab 2/4 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	62kDa
GeneID?Human?	GAB2/GAB4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GAB4; GRB2-associated-binding protein 4; GRB2-associated binder 2-like; GAB2-like; GRB2-associated binder 4; GRB2-associated-binding protein 2-like; Growth factor receptor bound protein 2-associated p
Background	similarity:Belongs to the GAB family.,similarity:Contains 1 PH domain.,