

FAM3D Polyclonal Antibody

Cat No: HR1AP5871

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | FAM3D Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,IHC-p,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 25kDa |
| GeneID?Human? | FAM3D |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | FAM3D; Protein FAM3D |
| Background | similarity:Belongs to the FAM3 family.,tissue specificity:Abundantly expressed in placenta and weakly expressed in small intestine., |