

EDEM2 Polyclonal Antibody

Cat No: HR1AP7866

For research use only

Overview

Product Name	EDEM2 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	EDEM2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	EDEM2; C20orf31; C20orf49; ER degradation-enhancing alpha-mannosidase-like 2
Background	ER degradation enhancing alpha-mannosidase like protein 2(EDEM2) Homo sapiens In the endoplasmic reticulum (ER), misfolded proteins are retrotranslocated to the cytosol and degraded by the proteasome in a process known as ER-associated degradation (ERAD). EDEM2 belongs to a family of proteins involved in ERAD of glycoproteins (Mast et al., 2005 [PubMed 15537790]).[supplied by OMIM, Mar 2008],