

DnaJC17 Polyclonal Antibody

Cat No: HR1AP7825

For research use only

Overview

Product Name	DnaJC17 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GeneID?Human?	DNAJC17
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	DNAJC17; DnaJ homolog subfamily C member 17
Background	similarity:Contains 1 J domain.,similarity:Contains 1 RRM (RNA recognition motif) domain.,