

CYTL1 Polyclonal Antibody

Cat No: HR1AP7761

For research use only

Overview

Product Name	CYTL1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	16kDa
GeneID?Human?	CYTL1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CYTL1; C4orf4; Cytokine-like protein 1; Protein C17
Background	cytokine like 1(CYTL1) Homo sapiens C17 is a cytokine-like protein specifically expressed in bone marrow and cord blood mononuclear cells that bear the CD34 (MIM 142230) surface marker (Liu et al., 2000 [PubMed 10857752]).[supplied by OMIM, Mar 2008].