

CPZ Polyclonal Antibody

Cat No: HR1AP8921

For research use only

Overview

Product Name	CPZ Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	73kDa
GeneID?Human?	CPZ
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CPZ; Carboxypeptidase Z; CPZ
Background	<p>carboxypeptidase Z(CPZ) Homo sapiens This gene encodes a member of the metalloproteinase family. This enzyme displays carboxypeptidase activity towards substrates with basic C-terminal residues. It is most active at neutral pH and is inhibited by active site-directed inhibitors of metalloproteinases. Alternative splicing in the coding region results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008].</p>