

CLN6 Polyclonal Antibody

Cat No: HR1AP7798

For research use only

Overview

Product Name	CLN6 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	CLN6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CLN6; Ceroid-lipofuscinosis neuronal protein 6; Protein CLN6
Background	ceroid-lipofuscinosis, neuronal 6, late infantile, variant(CLN6) Homo sapiens This gene is one of eight which have been associated with neuronal ceroid lipofuscinoses (NCL). Also referred to as Batten disease, NCL comprises a class of autosomal recessive, neurodegenerative disorders affecting children. The genes responsible likely encode proteins involved in the degradation of post-translationally modified proteins in lysosomes. The primary defect in NCL disorders is thought to be associated with lysosomal storage function. [provided by RefSeq, Oct 2008].