

CPXM2 Polyclonal Antibody

Cat No: HR1AP5785

For research use only

Overview

Product Name	CPXM2 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	CPXM2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CPXM2; Inactive carboxypeptidase-like protein X2
Background	caution:As it has lost active sites residues and zinc-binding sites it is unlikely to be catalytically active.,function:May be involved in cell-cell interactions.,similarity:Belongs to the peptidase M14 family.,similarity:Contains 1 F5/8 type C domain.,