

CHST1 Polyclonal Antibody

Cat No: HR1AP8922

For research use only

Overview

Product Name	CHST1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	49kDa
GeneID?Human?	CHST1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CHST1; Carbohydrate sulfotransferase 1; Galactose/N-acetylglucosamine/N-acetylglucosamine 6-O-sulfotransferase 1; GST-1; Keratan sulfate Gal-6 sulfotransferase; KS6ST; KSGal6ST; KSST
Background	carbohydrate sulfotransferase 1(CHST1) Homo sapiens This locus encodes a member of the keratin sulfotransferase family of proteins. The encoded enzyme catalyzes the sulfation of the proteoglycan keratin. [provided by RefSeq, Aug 2011],