

CEP76 Polyclonal Antibody

Cat No: HR1AP8714

For research use only

Overview

Product Name	CEP76 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	74kDa
GeneID?Human?	CEP76
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CEP76; C18orf9; Centrosomal protein of 76 kDa; Cep76
Background	centrosomal protein 76(CEP76) Homo sapiens This gene encodes a centrosomal protein which regulates centriole amplification by limiting centriole duplication to once per cell cycle. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2012],