

CAT-4 Polyclonal Antibody

Cat No: HR1AP8342

For research use only

Overview

Product Name	CAT-4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	70kDa
GeneID?Human?	SLC7A4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SLC7A4; Cationic amino acid transporter 4; CAT-4; CAT4; Solute carrier family 7 member 4
Background	function:Involved in the transport of the cationic amino acids (arginine, lysine and ornithine).,similarity:Belongs to the amino acid-polyamine-organocation (APC) superfamily. Cationic amino acid transporter (CAT) (TC 2.A.3.3) family.,