

BAI-2 Polyclonal Antibody

Cat No: HR1AP8035

For research use only

Overview

Product Name	BAI-2 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	BAI2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BAI2; Brain-specific angiogenesis inhibitor 2
Background	<p>adhesion G protein-coupled receptor B2(ADGRB2) Homo sapiens This gene encodes a a seven-span transmembrane protein that is thought to be a member of the secretin receptor family. The encoded protein is a brain-specific inhibitor of angiogenesis. The mature peptide may be further cleaved into additional products (PMID:20367554). Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2014],</p>