

## Ataxin-7L1 Polyclonal Antibody

Cat No: HR1AP6346

For research use only

### Overview

Product Name	Ataxin-7L1 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	ATXN7L1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ATXN7L1; ATXN7L4; KIAA1218; Ataxin-7-like protein 1; Ataxin-7-like protein 4
Background	similarity:Contains 1 SCA7 domain.,