

ARRDC3 Polyclonal Antibody

Cat No: HR1AP8031

For research use only

Overview

Product Name	ARRDC3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	46kDa
GeneID?Human?	ARRDC3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ARRDC3; KIAA1376; Arrestin domain-containing protein 3; TBP-2-like inducible membrane protein; TLIMP
Background	arrestin domain containing 3(ARRDC3) Homo sapiens This gene encodes a member of the arrestin family of proteins, which regulate G protein-mediated signaling. The encoded protein is thought to act as a regulator of breast cancer growth and progression by binding to a phosphorylated form of integrin beta4, a tumor-related antigen, targeting the integrin for internalization and degradation. [provided by RefSeq, Jul 2016],