

## Arrdc4 Polyclonal Antibody

Cat No: HR1AP9096

For research use only

### Overview

Product Name	Arrdc4 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	ARRDC4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ARRDC4; Arrestin domain-containing protein 4
Background	similarity:Belongs to the arrestin family.,