

ANKRD26 Polyclonal Antibody

Cat No: HR1AP6371

For research use only

Overview

Product Name	ANKRD26 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	ANKRD26
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ANKRD26; KIAA1074; Ankyrin repeat domain-containing protein 26
Background	ankyrin repeat domain 26(ANKRD26) Homo sapiens This gene encodes a protein containing N-terminal ankyrin repeats which function in protein-protein interactions. Mutations in this gene are associated with autosomal dominant thrombocytopenia-2. Pseudogenes of this gene are found on chromosome 7, 10, 13 and 16. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2011],