

Anoctamin-9 Polyclonal Antibody

Cat No: HR1AP6920

For research use only

Overview

Product Name	Anoctamin-9 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	ANO9
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ANO9; PIG5; TMEM16J; TP53I5; Anoctamin-9; Transmembrane protein 16J; Tumor protein p53-inducible protein 5; p53-induced gene 5 protein
Background	function:May act as a calcium-activated chloride channel.,sequence caution:Intron retention.,similarity:Belongs to the anoctamin family.,