

ANKRD20 Polyclonal Antibody

Cat No: HR1AP7383

For research use only

Overview

Product Name	ANKRD20 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	95kDa
GeneID?Human?	ANKRD20A3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ANKRD20A3; Ankyrin repeat domain-containing protein 20A3
Background	similarity:Contains 5 ANK repeats.,tissue specificity:Widely expressed; not found in kidney, nor in liver.,