

Angiotensin-L1 Polyclonal Antibody

Cat No: HR1AP6036

For research use only

Overview

Product Name	Angiotensin-L1 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	AMOTL1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	AMOTL1; Angiotensin-like protein 1
Background	<p>angiotensin like 1(AMOTL1) Homo sapiens The protein encoded by this gene is a peripheral membrane protein that is a component of tight junctions or TJs. TJs form an apical junctional structure and act to control paracellular permeability and maintain cell polarity. This protein is related to angiotensin, an angiotensin binding protein that regulates endothelial cell migration and capillary formation. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2014],</p>