

TFPI-2 Polyclonal Antibody

Cat No: HR1AP5000

For research use only

Overview

Product Name	TFPI-2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	13kDa
GeneID?Human?	TFPI2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TFPI2; Tissue factor pathway inhibitor 2; TFPI-2; Placental protein 5; PP5
Background	tissue factor pathway inhibitor 2(TFPI2) Homo sapiens This gene encodes a member of the Kunitz-type serine proteinase inhibitor family. The protein can inhibit a variety of serine proteases including factor VIIa/tissue factor, factor Xa, plasmin, trypsin, chymotrypsin and plasma kallikrein. This gene has been identified as a tumor suppressor gene in several types of cancer. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2012],