

FAM3D Polyclonal Antibody

Cat No: HR1AP9607

For research use only

Overview

Product Name	FAM3D Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	FAM3D
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Protein FAM3D
Background	similarity:Belongs to the FAM3 family.,tissue specificity:Abundantly expressed in placenta and weakly expressed in small intestine.,