

## Elmo1 Polyclonal Antibody

Cat No: HR1AP4871

For research use only

### Overview

Product Name	Elmo1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	83kDa
GeneID?Human?	ELMO1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ELMO1; KIAA0281; Engulfment and cell motility protein 1; Protein ced-12 homolog
Background	<p>engulfment and cell motility 1(ELMO1) Homo sapiens This gene encodes a member of the engulfment and cell motility protein family. These proteins interact with dedicator of cytokinesis proteins to promote phagocytosis and cell migration. Increased expression of this gene and dedicator of cytokinesis 1 may promote glioma cell invasion, and single nucleotide polymorphisms in this gene may be associated with diabetic nephropathy. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013],</p>