

ZnT-9 Polyclonal Antibody

Cat No: HR1AP5591

For research use only

Overview

Product Name	ZnT-9 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GeneID?Human?	SLC30A9
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>SLC30A9 (Solute Carrier Family 30 Member 9) is a Protein Coding gene. Diseases associated with SLC30A9 include Birk-Landau-Perez Syndrome. Among its related pathways are Metal ion SLC transporters. Gene Ontology (GO) annotations related to this gene include DNA binding transcription factor activity and ligand-dependent nuclear receptor transcription coactivator activity. Acts as a zinc transporter involved in intracellular zinc homeostasis (PubMed: 28334855). Functions as a secondary coactivator for nuclear receptors by cooperating with p160 coactivators subtypes. Plays a role in transcriptional activation of Wnt-responsive genes.</p>