

## ?-1,3-Gal-TL Polyclonal Antibody

Cat No: HR1AP4860

For research use only

### Overview

Product Name	?-1,3-Gal-TL Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	57kDa
GeneID?Human?	B3GALTL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	B3GALTL; B3GTL; Beta-1; 3-glycosyltransferase; Beta3Glc-T; Beta-3-glycosyltransferase-like
Background	beta 3-glycosyltransferase(B3GLCT) Homo sapiens The protein encoded by this gene is a beta-1,3-glycosyltransferase that transfers glucose to O-linked fucosylglycans on thrombospondin type-1 repeats (TSRs) of several proteins. The encoded protein is a type II membrane protein. Defects in this gene are a cause of Peters-plus syndrome (PPS).[provided by RefSeq, Mar 2009],