

Tollip Polyclonal Antibody

Cat No: HR1AP4736

For research use only

Overview

Product Name	Tollip Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	TOLLIP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TOLLIP; Toll-interacting protein
Background	toll interacting protein(TOLLIP) Homo sapiens This gene encodes a ubiquitin-binding protein that interacts with several Toll-like receptor (TLR) signaling cascade components. The encoded protein regulates inflammatory signaling and is involved in interleukin-1 receptor trafficking and in the turnover of IL1R-associated kinase. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2016],