

TIGD3 Polyclonal Antibody

Cat No: HR1AP4719

For research use only

Overview

Product Name	TIGD3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	TIGD3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TIGD3; Tigger transposable element-derived protein 3
Background	tigger transposable element derived 3(TIGD3) Homo sapiens The protein encoded by this gene belongs to the tigger subfamily of the pogo superfamily of DNA-mediated transposons in humans. These proteins are related to DNA transposons found in fungi and nematodes, and more distantly to the Tc1 and mariner transposases. They are also very similar to the major mammalian centromere protein B. The exact function of this gene is not known. [provided by RefSeq, Jul 2008],