

T3JAM Polyclonal Antibody

Cat No: HR1AP4658

For research use only

Overview

Product Name	T3JAM Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	64kDa
GeneID?Human?	TRAF3IP3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TRAF3IP3; T3JAM; TRAF3-interacting JNK-activating modulator; TRAF3-interacting protein 3
Background	TRAF3 interacting protein 3(TRAF3IP3) Homo sapiens The gene encodes a protein that mediates cell growth by modulating the c-Jun N-terminal kinase signal transduction pathway. The encoded protein may also interact with a large multi-protein assembly containing the phosphatase 2A catalytic subunit. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2013],