

## Sideroflexin-4 Polyclonal Antibody

Cat No: HR1AP4555

For research use only

### Overview

Product Name	Sideroflexin-4 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	SFXN4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SFXN4; BCRM1; Sideroflexin-4; Breast cancer resistance marker 1
Background	sideroflexin 4(SFXN4) Homo sapiens This gene encodes a member of the sideroflexin family. The encoded protein is a transmembrane protein of the inner mitochondrial membrane, and is required for mitochondrial respiratory homeostasis and erythropoiesis. Mutations in this gene are associated with mitochondriopathy and macrocytic anemia. Alternatively spliced transcript variants have been found in this gene. [provided by RefSeq, Jan 2014],