

Septin 14 Polyclonal Antibody

Cat No: HR1AP6951

For research use only

Overview

Product Name	Septin 14 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GeneID?Human?	Sep-14
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SEPT14; Septin-14
Background	septin 14(SEPT14) Homo sapiens SEPT14 is a member of the highly conserved septin family of GTP-binding cytoskeletal proteins implicated in membrane transport, apoptosis, cell polarity, cell cycle regulation, cytokinesis, and other cellular functions (Peterson et al., 2007 [PubMed 17922164]).[supplied by OMIM, Jul 2008],