

RIN3 Polyclonal Antibody

Cat No: HR1AP4492

For research use only

Overview

Product Name	RIN3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	107kDa
GeneID?Human?	RIN3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RIN3; Ras and Rab interactor 3; Ras interaction/interference protein 3
Background	Ras and Rab interactor 3(RIN3) Homo sapiens Summary: This protein encoded by this gene is a member of the RIN family of Ras interaction-interference proteins, which are binding partners to the RAB5 small GTPases. The protein functions as a guanine nucleotide exchange for RAB5B and RAB31. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2016],