

Ribosomal Protein S27L Polyclonal Antibody

Cat No: HR1AP4483

For research use only

Overview

Product Name	Ribosomal Protein S27L Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	18kDa
GeneID?Human?	RPS27L
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RPS27L; 40S ribosomal protein S27-like
Background	ribosomal protein S27 like(RPS27L) Homo sapiens This gene encodes a protein sharing 96% amino acid similarity with ribosomal protein S27, which suggests the encoded protein may be a component of the 40S ribosomal subunit. [provided by RefSeq, Jul 2008],