

Ribosomal Protein L26L Polyclonal Antibody

Cat No: HR1AP4471

For research use only

Overview

Product Name	Ribosomal Protein L26L Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	20kDa
GeneID?Human?	RPL26L1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RPL26L1; RPL26P1; 60S ribosomal protein L26-like 1
Background	ribosomal protein L26 like 1(RPL26L1) Homo sapiens This gene encodes a protein that shares high sequence similarity with ribosomal protein L26. Alternative splicing results in multiple transcript variants encoding the same protein. [provided by RefSeq, Dec 2015],