

# RGS12 Polyclonal Antibody

Cat No: HR1AP8154

For research use only

## Overview

Product Name	RGS12 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	160kDa
GeneID?Human?	RGS12
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RGS12; Regulator of G-protein signaling 12; RGS12
Background	<p>regulator of G-protein signaling 12(RGS12) Homo sapiens This gene encodes a member of the 'regulator of G protein signaling' (RGS) gene family. The encoded protein may function as a guanosine triphosphatase (GTPase)-activating protein as well as a transcriptional repressor. This protein may play a role in tumorigenesis. Multiple transcript variants encoding distinct isoforms have been identified for this gene. Other alternative splice variants have been described but their biological nature has not been determined. [provided by RefSeq, Jul 2008].</p>