

Recoverin Polyclonal Antibody

Cat No: HR1AP8129

For research use only

Overview

Product Name	Recoverin Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	23kDa
GeneID?Human?	RCVRN
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RCVRN; RCV1; Recoverin; Cancer-associated retinopathy protein; Protein CAR
Background	<p>recoverin(RCVRN) Homo sapiens This gene encodes a member of the recoverin family of neuronal calcium sensors. The encoded protein contains three calcium-binding EF-hand domains and may prolong the termination of the phototransduction cascade in the retina by blocking the phosphorylation of photo-activated rhodopsin. Recoverin may be the antigen responsible for cancer-associated retinopathy. [provided by RefSeq, Jul 2008],</p>