

PRSS33 Polyclonal Antibody

Cat No: HR1AP4376

For research use only

Overview

Product Name	PRSS33 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	26kDa
GeneID?Human?	PRSS33
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PRSS33; Serine protease 33; Serine protease EOS
Background	<p>function:Serine protease that has amidolytic activity, cleaving its substrates before Arg residues.,induction:Up-regulated by phorbol ester PMA.,similarity:Belongs to the peptidase S1 family.,similarity:Contains 1 peptidase S1 domain.,tissue specificity:Predominantly expressed in macrophages. Present in the spleen, small and large intestine, lung and brain (at protein level). Highly expressed in peripheral leukocytes, ovary, retina, spleen and stomach. Moderately expressed in thymus, uterus and platelets, as well as some brain tissues, such as thalamus and fetal brain.,</p>