

RFPL4A Polyclonal Antibody

Cat No: HR1AP4455

For research use only

Overview

Product Name	RFPL4A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	RFPL4A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RFPL4A; RFPL4; RNF210; Ret finger protein-like 4A; RING finger protein 210
Background	The ret finger protein-like 4 gene, Rfpl4, encodes a putative E3 ubiquitin-protein ligase expressed in adult germ cells.