



PWWP2B Polyclonal Antibody

Cat No: HR1AP6109

For research use only

Overview

Product Name	PWWP2B Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	64kDa
GeneID?Human?	PWWP2B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PWWP2B; PWWP2; PWWP domain-containing protein 2B
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,sequence caution:Translation N-terminally shortened.,similarity:Contains 1 PWWP domain.,