

Photomedin-2 Polyclonal Antibody

Cat No: HR1AP4306

For research use only

Overview

Product Name	Photomedin-2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	84kDa
GeneID?Human?	OLFML2B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	OLFML2B; Olfactomedin-like protein 2B; Photomedin-2
Background	PTM:O-glycosylated and N-glycosylated.,similarity:Contains 1 olfactomedin-like domain.,subunit:Homodimer. Binds to heparin and chondroitin sulfate E.,