

Photomedin-1 Polyclonal Antibody

Cat No: HR1AP6107

For research use only

Overview

Product Name	Photomedin-1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	73kDa
GeneID?Human?	OLFML2A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	OLFML2A; Olfactomedin-like protein 2A; Photomedin-1
Background	PTM:May be cleaved at Lys-295 after secretion.,PTM:O-glycosylated but not N-glycosylated.,similarity:Contains 1 olfactomedin-like domain.,subunit:Homodimer. Binds to heparin and chondroitin sulfate E.,