

NARG1 Polyclonal Antibody

Cat No: HR1AP4009

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | NARG1 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,IHC-p,IF,ELISA |
| Species Reactivity | Human,Mouse |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 100kDa |
| GeneID?Human? | NAA15 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | NAA15; GA19; NARG1; NATH; TBDN100; N-alpha-acetyltransferase 15; NatA auxiliary subunit; Gastric cancer antigen Ga19; N-terminal acetyltransferase; NMDA receptor-regulated protein 1; Protein tubedown |
| Background | N(alpha)-acetyltransferase 15, NatA auxiliary subunit(NAA15) Homo sapiens This gene encodes a protein of unknown function. However, similarity to proteins in yeast and other species suggests that this protein may be an N-acetyltransferase. [provided by RefSeq, Jul 2008], |