

MRP-L9 Polyclonal Antibody

Cat No: HR1AP3956

For research use only

Overview

Product Name	MRP-L9 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	MRPL9
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MRPL9; 39S ribosomal protein L9; mitochondrial; L9mt; MRP-L9
Background	mitochondrial ribosomal protein L9(MRPL9) Homo sapiens This is a nuclear gene encoding a protein component of the 39S subunit of the mitochondrial ribosome. Alternative splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 8. [provided by RefSeq, Jul 2014],