

# MRP-L11 Polyclonal Antibody

Cat No: HR1AP3938

For research use only

## Overview

Product Name	MRP-L11 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	21kDa
GeneID?Human?	MRPL11
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MRPL11; CGI-113; 39S ribosomal protein L11; mitochondrial; L11mt; MRP-L11
Background	mitochondrial ribosomal protein L11(MRPL11) Homo sapiens This nuclear gene encodes a 39S subunit component of the mitochondrial ribosome. Alternative splicing results in multiple transcript variants. Pseudogenes for this gene are found on chromosomes 5 and 12. [provided by RefSeq, May 2014],